



COVID-19 Ab test

Anti-viral antibody IgM/IgG test

- Easy to use
- **Fast and precise result (10 mins)**
- Small POC device
- Room temperature storage
- Anti-viral IgG/IgM test
- Automated (AFIAS) or manual(ichroma II) platform

3 simple steps



1 Collect capillary blood






2 Place the C-tip™ in the cartridge tip hole & run the test



3 Press START to run the test.



High sensitivity (Time resolved fluorescence) immunoassay



Product name	AFIAS COVID-19 Ab	ichroma COVID-19 Ab
Assay type	TRFLFA*, IgM/IgG	
Sample type	Whole blood/ Serum/plasma	
Sample prep/ extraction	Fingertip blood(C-tip)/ tube blood	
Target	Anti-viral IgM/IgG	
Platform	AFIAS-1/6 (Automated)	ichroma II (Manual)
Unit	24 kit/box	25 kit/box